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Operating instructions

STEAM

for Dentsply Sirona chemical indicators

REF	Product Code	Application	Norm EN ISO 11140-1	Sterilization process	Stated Value	Content [pieces]
6638220	DS chemical indicator class 5 (500x)	Self-adhesive data labels with indicator to monitor sterilization processes and subsequent patient-related documentation	Type 5	Steam	134 °C, 3 min	500

Application

All Dentsply Sirona chemical indicators contain an integrating indicator that allows reliable information about the quality of the sterilization at the location where it was placed. Temperature and pressure versus time are easily monitored by the sterilizer. However, air or non condensable gases (NCG) in the sterilization chamber cannot easily be detected by those physical test methods. The presence of NCG negatively influences the result of the sterilization process. This may occur by insufficient air removal, leakages in the door seals or valves or in most cases NCG enter via the external sterilization gas supply into the sterilization chamber. Especially the combination of NCG, low temperature and dense packs may result in faulty sterilization processes.

If only cross contamination on surfaces of the instruments should be prevented, the instruments are sterilized unpackaged. In special combination autoclaves that clean, lubricate and sterilize dental instruments, the indicator is placed unpacked into the holder inside the chamber.

If the Dentsply Sirona chemical indicator is placed into the package/container, where the sterilization agent penetration is most difficult, the security is given that sterile conditions were reached when the indicator changed its colour according to the description on the indicator. If the colour is not changing according to the description on the indicator, this is an indication of presence of NCG in the pack resulting in too low concentration of the sterilant and/or an insufficient temperature-time window in the pack.

Packaging monitoring indicators should only be used if solid and porous goods are sterilized. They are not recommended if hollow devices are sterilized. In this case they may give false-positive information and Batch Monitoring Systems (BMS) based on an adequate hollow test Process Challenge Device (PCD) have to be used.

Unfortunately, the result of the sterilization process can only be determined after opening the pack. Therefore, it is important to get the information on a successful sterilization process shortly after the process without opening the packs or containers. The batch monitoring PCD simulates the worst case penetration conditions in a pack. The integrating indicator placed inside the PCD can be taken out immediately after the sterilization process. The operator has the information about the efficiency of the sterilization process. Faulty batches are immediately detected and do not get into the operation room.

Handling Information

1. Select the correct indicator for the sterilization process (see table on page 1).
2. Place the indicator in the middle of the package or container (most difficult location to sterilize). Alternatively place the indicator unpacked inside the chamber. In DAC Universal combination autoclaves the indicator is fixed inside the indicator holder.
3. Run the sterilization program. At the end of the sterilization cycle or after opening of the package, remove the indicator and check the result:
 - If the indicator or all indicator bars have changed its colour to the final colour (see label for colour change) the sterilization was successful.
 - If the colour did not change to the final colour, the success of the sterilization process cannot be guaranteed. The content of the package should not be used in the operation theatre.
4. All indicators are adhesive and can be documented onto the patient documentation sheet to be able to check the result later on in case a nosocomial infection occurs to exclude incomplete sterilization.

Storage and Disposal

1. For longer periods store all indicators in the original package.
2. Store indicators always between 5-30°C or 41-86°F and a humidity of 5 - 80% RH.
3. The vapour of chemicals especially hydrogen peroxide may change the indicator before or after sterilization. Therefore, do not store them together with other chemicals.
4. The indicators should not be used after expiry date. They may be disposed off with normal waste.

Safety Precautions

Monitoring indicators do not replace the basic validation according to national and/ or international standards. Follow the instructions of your local requirements.